



DNA Ligation Kit

Instruction Manual

Catalog #203003

Revision B

Research Use Only. Not for Use in Diagnostic Procedures.

203003-12



Agilent Technologies

LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Agilent. Agilent shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

ORDERING INFORMATION AND TECHNICAL SERVICES

Email

techservices@agilent.com

World Wide Web

www.genomics.agilent.com

Telephone

Location	Telephone
United States and Canada	800 227 9770
Austria	01 25125 6800
Benelux	02 404 92 22
Denmark	45 70 13 00 30
Finland	010 802 220
France	0810 446 446
Germany	0800 603 1000
Italy	800 012575
Netherlands	020 547 2600
Spain	901 11 68 90
Sweden	08 506 4 8960
Switzerland	0848 8035 60
UK/Ireland	0845 712 5292
All Other Countries	Please visit www.genomics.agilent.com and click Contact Us

DNA Ligation Kit

CONTENTS

Materials Provided	1
Storage Conditions	1
Introduction	1
Preprotocol Considerations	2
Protocol	2
Ligation of Insert DNA into Lambda Vector Arms	2
Ligation of Insert DNA into Plasmid Vectors	3
Preparation of Media and Reagents	4
MSDS Information	4

DNA Ligation Kit

MATERIALS PROVIDED

Note *Aliquot the 10 mM rATP into smaller volumes after initial thawing to avoid multiple freeze–thaw cycles.*

Materials provided	Quantity	Storage temperature
10 mM rATP (pH 7.5) in sterile water	4 × 250 µl (1 ml total volume)	–20°C ^a
cl857 wild-type lambda control DNA, <i>Hind</i> III digested	10 µg	–20°C
pUC18 plasmid control DNA, <i>Bam</i> H I digested	10 µg	–20°C
T4 DNA ligase (4 U/µl)	300 U	–20°C
10× ligase buffer ^b	1 ml	–20°C

^a For long-term storage, store at –80°C.

^b See *Preparation of Media and Reagents*.

STORAGE CONDITIONS

10 mM rATP (pH 7.5): –20°C (–80°C for long-term storage)

cl857 Wild-Type Lambda Control DNA: –20°C

pUC18 Plasmid Control DNA: –20°C

T4 DNA Ligase: –20°C

10× Ligase Buffer: –20°C

INTRODUCTION

The DNA Ligation Kit contains the reagents necessary to perform both sticky and blunt-end ligations. These reagents are guaranteed to be DNase free while providing optimal ligation efficiency. The kit contains reagents for 150 ligations and is optimized for use with Agilent's lambda vector arms and plasmid vectors.

PREPROTOCOL CONSIDERATIONS

- ♦ When setting up ligations, **gently** mix the ligation reaction with the pipet tip after adding the ligase.
- ♦ T4 DNA ligase is not inhibited by transfer RNA (tRNA), but is strongly inhibited by NaCl concentrations >150 mM. To remove NaCl, wash DNA pellets with 80% (v/v) alcohol prior to resuspension.
- ♦ The efficiency of any ligation reaction depends on the integrity of the cohesive ends being ligated and on the quality of the vector DNA.
- ♦ Plasmid blunt-end ligations may be enhanced by reducing the final rATP concentration by half and lengthening the ligation time.
- ♦ When performing any ligation experiment, it is important to include a control ligation of the prepared vector alone. This will provide a good evaluation of the background associated with the vector in the experiment.

The following formula allows easy calculation of picomoles per end:

$$\frac{2 \times 10^6}{660 \times \text{number of base pairs}} = \text{picomole ends per microgram of DNA}$$

PROTOCOL

Ligation of Insert DNA into Lambda Vector Arms

Control Ligation

To verify the activity of the components, perform the following control ligation using 1 µl of the *cI857* wild-type lambda control DNA.

1. Add the following components to a microcentrifuge tube:

1 µg (1–2 µl) of *cI857* wild-type lambda control DNA
0.5 µl of 10× ligase buffer
0.5 µl of 10 mM rATP (pH 7.5)
X µl of sterile water for a final volume of 4.5 µl

Then add

0.5 µl of T4 DNA ligase (4 U/µl)

Note *Keep the final glycerol content of the T4 DNA ligase below 5%.*

2. Incubate the control ligation reaction at 12°C for 4–6 hours or at 4°C overnight.

3. Prepare a 0.6% (w/v) agarose gel and load 1 μl of unligated *cI857* wild-type lambda control DNA and 5 μl of the control ligation reaction. After running the gel, compare the unligated lanes with the ligated lanes to verify evidence of activity.

Sample Ligation

1. Add the following components to a microcentrifuge tube:

1 μg of prepared lambda vector arms
 X μg of insert DNA^{||}
0.5 μl of 10 \times ligase buffer
0.5 μl of 10 mM rATP (pH 7.5)
 X μl of sterile water for a final volume of 4.5 μl

Then add

0.5 μl of T4 DNA ligase (4 U/ μl)

Note *Keep the final glycerol content of the T4 DNA ligase below 5%.*

2. Incubate the sample ligation reaction at 12°C for 4–6 hours or at 4°C overnight.

^{||} When ligating into lambda vector arms, we recommend using an equimolar insert-to-vector ratio in order to prevent multiple inserts.

Ligation of Insert DNA into Plasmid Vectors

Control Ligation

To verify the activity of the components, perform the following control ligation using 1 μl of pUC18 plasmid control DNA.

1. Add the following components to a microcentrifuge tube:

1 μg (1–2 μl) of pUC18 plasmid control DNA
1 μl of 10 \times ligase buffer
0.5 μl of 10 mM rATP (pH 7.5)
 X μl of sterile water for a final volume of 9.5 μl

Then add

0.5 μl of T4 DNA ligase (4 U/ μl)

Note *Keep the final glycerol content of the T4 DNA ligase below 5%.*

2. Incubate the control ligation reaction at 12°C for 4–6 hours or at 4°C overnight.

3. Prepare a 0.6% (w/v) agarose gel and load 1 μ l of unligated pUC18 plasmid control DNA and 10 μ l of the control ligation reaction. After running the gel, compare the unligated lanes with the ligated lanes to verify evidence of activity.

Sample Ligation

1. Add the following components to a microcentrifuge tube:

X μ g of vector DNA
 X μ g of insert DNA[¶]
1 μ l of 10 \times ligase buffer
1 μ l of 10 mM rATP (pH 7.5)
0.5 μ l of T4 DNA ligase (4 U/ μ l)
 X μ l of sterile water for a final volume of 10 μ l

Note *Keep the final glycerol content of the T4 DNA ligase below 5%.*

2. Incubate the sample ligation reaction at 12°C for 4–6 hours or at 4°C overnight.

[¶] Each vector–insert combination requires unique ligation conditions. To obtain optimum results, we recommend testing various vector-to-insert ratios to ensure the highest ligation efficiency (e.g., insert-to-vector ratios from 2:1 to 10:1).

PREPARATION OF MEDIA AND REAGENTS

10 \times Ligase Buffer

500 mM Tris-HCl (pH 7.5)
70 mM MgCl₂
10 mM dithiothreitol (DTT)

Note *rATP is added separately in the ligation reaction*

MSDS INFORMATION

Material Safety Data Sheets (MSDSs) are provided online at <http://www.genomics.agilent.com>. MSDS documents are not included with product shipments.