# SURESELECT<sup>XT HS</sup> TARGET ENRICHMENT 超微量目標外顯子組富集系統

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
Figure 1. SureSelect<sup>XT HS</sup> optimized performance for FFPE with the majority of bases covered at 20X across decreasing sample quality. % bases covered by at least 20X read depth with increasing sequencing depth (4-8 million - 2x100 bp reads) is plotted against gDNA quality (qPCR-based gDNA quality (ddCq) and is calculated using Agilent's FFPE QC Assay kit).


Figure 2. SureSelect<sup>XT HS</sup> outperforms other library prep solutions at 100X coverage, especially as sample quality decreases: % bases covered by at least 100X read depth is plotted against gDNA of decreasing quality (qPCR-based gDNA quality (ddCq) and is calculated using Agilent's FFPE QC Assay kit (based on 1000X sequencing depth with 8 million - 2x100bp reads).

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
**Figure 3. Detect down to 1% allele frequency.**  
correlation of observed to expected frequency down to as little as 1% VAF (HD200 Quantitative Multiplex Reference Standard and ClearSeq Comprehensive Cancer Panel).

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
### Basic molecular barcode (MBC) analysis


- Step ① Align reads
- Step ② Group read pairs to designed probes based on read start-stop position
- Step ③ Group reads with an identical molecular barcode sequence for each probe
- Step ④ Consolidate read information to one read per molecule (remove PCR duplicates)

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
### Standard deduplication

◆ SNVs verified by Horizon ◆ Other non-reference bases detected but not on the list provided by Horizon

### MBC deduplication

**Figure 4. Improved data quality with molecular barcodes (MBC) analysis.** Targeted sequencing of SureSelect<sup>XT HS</sup> libraries from Horizon cfDNA reference and a custom 164 kb SureSelect panel achieved a median coverage of 808(standard) / 906(MBC) after deduplication. 216 non-reference alleles were detected at a variant allele frequency of less than 2% with standard deduplication compared to 18 using the MBC. This represents a 92% reduction in false calls. Courtesy of Dr. L.J. Barber and Dr. M. Gerlinger, Centre for Evolution and Cancer, The Institute of Cancer Research, London, UK.

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### Create a Targeted Capture Panel

#### Genomic locations of interest

Design target sequences in SureDesign

#### Capture Library (SureSelect or ClearSeq)

### Library Preparation

#### Genomic DNA samples (10 ng - 200 ng DNA per sample)

Shear DNA

#### DNA fragments of 150 - 200 bp

Prepare samples using SureSelect<sup>XT HS</sup> Library Prep Kit

#### Molecular barcode & indexed DNA libraries

PCR amplify

#### Prepared DNA libraries

Hybridize using SureSelect<sup>XT HS</sup> reagents and protocol

#### Capture Library/prepared DNA library hybrids


Capture hybrids on magnetic beads

#### Captured libraries

PCR amplify

Pool captured samples (Optional)

#### SureSelect-enriched indexed NGS samples

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SureSelect Library Prep Solutions for Cancer and Constitutional Applications					
Product Name	SureSelect <sup>XT HS</sup>	SureSelect <sup>XT Low Input</sup>	SureSelect <sup>XT</sup>	SureSelect <sup>XT2</sup>	SureSelect <sup>XT</sup>
DNA Input	10 ng - 200 ng	10 ng - 200 ng	200 ng - 3 µg	100 ng - 1 µg	50 ng
Turnaround Time	8 hr	8 hr	1.5 day	1.5 day	7 hours
Covaris Needed	Yes	Yes	Yes	Yes	No
Library Complexity	Highest	Highest	High	Medium	Medium
Unique Features	FFPE optimized Molecular barcodes Mastermixed reagents Samples indexed prior to capture eliminating concern of cross sample contamination	FFPE optimized Molecular barcodes (optional) Mastermixed reagents Samples indexed prior to capture eliminating concern of cross sample contamination	Compatible with FFPE samples Robust variant identification	Pre-capture pooling Mastermix reagents	Transposase-based Mastermix reagents Whole genome sequencing and target enrichment compatible For intact DNA only
Key Benefits	High sensitivity for ≤1% VAF	192 sample indexes	High-complexity libraries for rare allele detection	Cost-effective	Covaris-free workflow

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© Agilent Technologies, Inc. 2017