



## TESTING COLONIES FOR ADENOVIRUS RECOMBINANTS

1. Examine the plates containing the pUC18 transformation to calculate the transformation efficiency (expect  $\geq 1 \times 10^8$  cfu/ $\mu$ g).
2. Compare the colonies on the linearized transfer vector alone (control) to the colonies on the recombination plates. The control plates should contain colonies of uniform size, arising from recircularized transfer vector. The transformants on the plates containing the adenoviral DNA recombinants will appear as two populations: normal-sized and tiny.
3. Pick 10 of the smallest, best-isolated colonies from the recombination plates (transfer vector plus gene of interest) into 2-ml cultures of LB broth containing the appropriate antibiotic (see *Preparation of Media and Reagents*).
4. Incubate at 37°C overnight while shaking at 225–250 rpm.
5. Prepare miniprep DNA from the overnight cultures using method of choice. The final volume of miniprep DNA should be 50  $\mu$ l and the DNA should be resuspended in sterile dH<sub>2</sub>O.

**Note** Do not store the BJ5183 cultures after overnight growth as undesired recombinants can be generated. Prepare plasmid miniprep DNA first thing in the morning.

6. Cut 5  $\mu$ l of the miniprep DNA with restriction enzymes that are diagnostic for the recombination event and run the digest on a 0.8% agarose TAE gel (see *Preparation of Media and Reagents*) next to 5  $\mu$ g of uncut plasmid. It is also recommended that 5  $\mu$ l of the miniprep DNA be cut with a restriction enzyme that will cleave somewhere within the gene of interest to confirm maintenance of the insert in the recombined adenovirus plasmid.

**Note** Reserve a small amount of each plasmid sample for amplifying by transformation in a subsequent step.

Once the construction of the recombinant adenovirus plasmid(s) has been confirmed, amplify the plasmid stock by transforming competent bacterial cells and purifying maxiprep DNA from these. It is not necessary to use BJ5183 cells for this as no recombination event takes place during amplification and many other bacterial strains give higher yields of DNA. Following amplification, packaged adenovirus can be produced by transfecting human cell lines with the linearized adenoviral DNA. Packaged adenovirus can then be used in gene expression studies.<sup>2</sup>

## PREPARATION OF MEDIA AND REAGENTS

LB Broth (per Liter)	LB Agar (per Liter)	LB-Ampicillin Agar (per Liter)	1 × TAE Buffer
10 g of NaCl 10 g of tryptone 5 g of yeast extract Add deionized H <sub>2</sub> O to a final volume of 1 liter Adjust to pH 7.0 with 5 N NaOH Autoclave Cool to 55°C Add antibiotic (if required)	10 g of NaCl 10 g of tryptone 5 g of yeast extract 20 g of agar Adjust pH to 7.0 with 5 N NaOH Add deionized H <sub>2</sub> O to a final volume of 1 liter Autoclave Cool to 55°C Add antibiotic (if required) Pour into petri dishes (~25 ml/100-mm dish)	1 liter of LB agar, autoclaved Cool to 55°C Add 10 ml of 10 mg/ml filter-sterilized ampicillin Pour into petri dishes (~25 ml/100-mm plate)	40 mM Tris-acetate 1 mM EDTA

## REFERENCES

1. Hanahan, D. (1983) *J. Mol. Biol.* 166: 557–580.
2. He, T., Zhou, S., DaCosta, L., Yu, J., Kinzler, K., and Vogtlestein, B. (1998) *Proc. Natl. Acad. Sci. USA* 95:2509–2514.

## QUALITY CONTROL TESTING

Following the protocol above, BJ5183 electroporation competent cells are transformed with and without pUC18 control plasmid DNA. One microliter of pUC18 DNA (10 pg/ $\mu$ l) is used to transform 40  $\mu$ l of BJ5183 cells. After electroporation, the cells are resuspended in 1 ml of LB broth and allowed to recover for 1 hour at 37°C with shaking. Following transformation, 5- $\mu$ l samples of the culture are plated in duplicate on LB agar plates with 100  $\mu$ g/ml of ampicillin. The plates are incubated overnight at 37°C. The efficiency is calculated based on the average number of colonies per plate. In addition, the BJ5183 electroporation competent cells are tested for adenovirus recombination proficiency using Stratagene's AdEasy® vector system.

## LIMITED PRODUCT WARRANTY

This warranty limits our liability to replacement of this product. No other warranties of any kind, express or implied, including without limitation, implied warranties of merchantability or fitness for a particular purpose, are provided by Stratagene. Stratagene shall have no liability for any direct, indirect, consequential, or incidental damages arising out of the use, the results of use, or the inability to use this product.

## ENDNOTES

Falcon® is a registered trademark of Beckton-Dickinson and Company.  
AdEasy is a registered trademark of Johns Hopkins University.

**Stratagene Technical Services**  
USA/Canada (Toll-free) 800 894 1304  
Europe (Toll-free) 00800 7400 7400  
**Email** tech\_services@stratagene.com  
**World Wide Web** www.stratagene.com